

MM - Measuring Microscope

Total Magnification = Optical Magnification x Digital Magnification

Optical Magnification = Objective Magnification x C-Mount Adapter Magnification

Digital Magnification = Screen Size(Diagonal of the **monitor** in mm) / Sensor Size(Diagonal of the **camera** in mm)

Optical Magnification

C-Mount Adapter

Eyepiece


Objective lens



Digital Magnification

Sensor(Camera)

Screen(Monitor)



Objective	C-mount	Monitor	Camera	Total
0.7	0.5	609.6	4.0	53.3
1.0	0.5	609.6	4.0	76.2
1.5	0.5	609.6	4.0	114.3
2.0	0.5	609.6	4.0	152.4
2.5	0.5	609.6	4.0	190.5

※ ex1) **Camera Size(1/4")** = SQRT(2.4²)+(3.2²)

※ ex2) **Monitor Size(24")** = 24×25.4mm

Sensor	Monitor Size						
	9"	13"	15"	17"	19"	23"	27"
1/4"	57.2	82.6	95.25	107.95	120.65	146.05	171.5
1/3"	38.1	55.0	63.5	72.0	80.43	97.37	114.3
1/2"	28.6	41.3	47.6	54.0	60.33	73.03	85.7
2/3"	20.8	30.0	34.6	39.3	43.87	53.11	62.3

※ VIDEO Magnification = Monitor diagonal length / CCD diagonal length


Camera Size(diagonal) Calculation

$$1/4" = \text{SQRT}(2.4^2 + (3.2^2)) = 4.0$$

$$1/3" = \text{SQRT}(3.6^2 + (4.8^2)) = 6.0$$

$$1/2" = \text{SQRT}(4.8^2 + (6.4^2)) = 8.0$$

$$2/3" = \text{SQRT}(6.6^2 + (8.8^2)) = 11.0$$



ex) Maximum magnification??

Objective Magnification = 4.5x

C-Mount Adapter Magnification = 50x


CCD Camera = 1/4"

Monitor = 27"

$$\rightarrow (4.5 \times 50 \times 685.8) / 4.0 = 38,576.3 \times$$

C-Mount Adapter

: The Nikon C-DA C-Mount Camera Adapter 1x is a C-Mount Adapter used to couple a digital camera to the intermediate tube of a Nikon microscope so the camera can properly capture the image projected from the microscope



IM - Industrial Microscope



IM - Stereo Microscope

